

ROTARY PENDULUM MODELING

Topics Covered

- Interact with the **QUANSER®** QUBE-Servo 2 Rotary Pendulum system.
- Configure sensor and actuator gains to match model conventions.

Prerequisites

- Hardware Interfacing laboratory experiment.

1 Background

The rotary pendulum system, also known as the Furuta Pendulum, is a classic system often used to teach modeling and control in physics and engineering. The convention used for the modeling of the QUBE-Servo 2 rotary pendulum is shown in Figure 1.1.




Figure 1.1: Model conventions of rotary pendulum

The rotary arm, which is attached to the motor pivot, is denoted by the variable θ and the pendulum angle, attached to the end of the rotary arm, is denoted by α . Note the following conventions:

- Angle α is defined as the *inverted pendulum angle*, i.e. the angle with respect to the upright vertical position where $\alpha = 0$ means the pendulum is perfectly upright. It is expressed mathematically using

$$\alpha = \alpha_{full} \bmod 2\pi - \pi. \quad (1.1)$$

where α_{full} is the pendulum angle measured by the encoder, i.e. the continuous angle measurement defined as zero when pendulum is in the downward configuration.

- Both angles are defined as positive when rotated in the counter-clockwise (CCW) direction.
- When a positive voltage is applied to the motor, the rotary arm moves in the positive CCW direction.

The goal is to design a system model that follows these conventions. The Hardware Interfacing laboratory experiment introduced the DC motor and encoders on the QUBE-Servo 2 system. The pendulum angle is also measured using an encoder.

2 In-Lab Exercises

In this lab, we will make a **SIMULINK®** model using **QUARC®** blocks to drive the dc motor and measure both the rotary arm and pendulum angles - similarly as shown in Figure 2.1.






Figure 2.1: **SIMULINK®** model used with **QUARC®** to drive motor and read angles on QUBE-Servo 2 Rotary Pendulum system

1. Using the **SIMULINK®** model you made in the Hardware Interfacing laboratory experiment, do the following:
 - Configure the HIL Read Encoder Timebase block to read from both channels #0 and #1. The pendulum is measured on channel #1.
 - Setup the encoder gains on each channel to read the angles in radians (instead of degrees as in the lab).
 - Connect the the measured angles to scopes, but display them in degrees. You can do this using the *Simulink Extras | Transformation* blocks or just by adding another Gain block that converts radians to degrees.
 - If desired, add a Slider Gain block between the Constant and Analog Output block. Make sure the Constant source block is set to 1 in this case.
2. Build and run the **QUARC®** controller.
3. Rotate the rotary arm and pendulum counter-clockwise and examine the response on the scopes. Example responses are shown in Figure 2.2. Do the measured angles follow the modeling conventions given in Section 1



(a) Rotary Arm



(b) Pendulum

Figure 2.2: Measured rotary arm and pendulum angles

4. Apply a small voltage (e.g. 0.5 V) to the motor by changing the Constant block connected to Analog Output. Does this adhere to the modeling conventions?
5. Modify the Simulink model such that the measured angles and applied voltage follow by the modeling conventions. Briefly list any changes made.
6. Add modulus and bias blocks, as shown in Figure 2.3, to measure *inverted pendulum angle*, defined as Equation 1.1.




Figure 2.3: **SIMULINK®** modulus and bias blocks

7. Make sure the pendulum is hanging in the downward position and lies motionless before starting the controller.
8. Build and run the **QUARC®** controller.
9. Rotate the pendulum to the upright vertical position and ensure the angle is measured correctly and it follows the model conventions shown in Figure 1.1. Capture the response of the pendulum being raised to the inverted position. Explain what the bias and modulus functions do?
10. Stop the **QUARC®** controller.

© 2020 Quanser Inc., All rights reserved.

Quanser Inc.
119 Spy Court
Markham, Ontario
L3R 5H6
Canada
info@quanser.com
Phone: 1-905-940-3575
Fax: 1-905-940-3576

Printed in Markham, Ontario.

For more information on the solutions Quanser Inc. offers, please visit the web site at:
<http://www.quanser.com>

This document and the software described in it are provided subject to a license agreement. Neither the software nor this document may be used or copied except as specified under the terms of that license agreement. Quanser Inc. grants the following rights: a) The right to reproduce the work, to incorporate the work into one or more collections, and to reproduce the work as incorporated in the collections, b) to create and reproduce adaptations provided reasonable steps are taken to clearly identify the changes that were made to the original work, c) to distribute and publically perform the work including as incorporated in collections, and d) to distribute and publicly perform adaptations. The above rights may be exercised in all media and formats whether now known or hereafter devised. These rights are granted subject to and limited by the following restrictions: a) You may not exercise any of the rights granted to You in above in any manner that is primarily intended for or directed toward commercial advantage or private monetary compensation, and b) You must keep intact all copyright notices for the Work and provide the name Quanser Inc. for attribution. These restrictions may not be waved without express prior written permission of Quanser Inc.