

Quick Installation Guide: QUARC™ on NI™ myRIO

STEP 1 Install MATLAB® and Required Add-Ons

QUARC™ supports 64-bit Microsoft® Windows®.

Ensure one of [supported MATLAB](#) versions is installed on the computer with the following required add-ons accompanying the corresponding MATLAB version:

- **Simulink®**
- **Simulink Coder™**
- **MATLAB Coder™** (required by Simulink Coder)
- **Control System Toolbox™** (Optional add-on, but highly recommended as used by most of Quanser's control laboratories)


STEP 2 Install NI CompactRIO™ Device Drivers

Install NI CompactRIO (version 18.0 or 19.6), which can be obtained from the [National Instruments™](#).

Select the default features, and ensure the NI Measurement & Automation Explorer (MAX) is selected to install on the host computer (NI MAX is required for configuring and installing QUARC on the myRIO).

STEP 3 Install QUARC on Windows


A



Uninstall any previous version of QUARC that may be present on the computer.

Do so by launching the *Programs and Features* dialog or the *Apps & features* dialog depending on which Windows version you have.

B



1. An internet connection is required during the QUARC installation process. Download the QUARC web installer executable using the link provided in the confirmation email that you received.
2. Run the QUARC installer (i.e., install_quarc.exe). The QUARC installation screen should appear.

The installer automatically checks if there is a newer update ready to download.

Tip: To find tips for each installation window, hover the mouse cursor




on the blue information icon on the lower left corner or



click the question mark icon on the upper right corner for more details from the installation guide.

C




If an update is available, choose to *Update* to the latest version [for free].

D



Read over the license agreement displayed in the Quanser License Agreement window.

E



Enter the location of the QUARC license file provided in the confirmation email and click *Next* to continue.

F



The installer will automatically scan the software environment on the host PC to ensure it meets the requirements for QUARC on the myRIO (e.g., MATLAB, NI MAX).


To start the installation immediately using the default settings click Install and skip to Step I. Otherwise click Options to customize the installation and continue.

G



Provide the destination folder where QUARC will be installed, and click *Next* to continue.

H



Choose the features to be installed. If you don't have the proper license required for a feature, the feature will be disabled and greyed out.

At minimum, make sure to select **myRIO Cross-Compilers** and **Simulink Development Environment**.

Note: If you don't have NI MAX installed, the *myRIO Cross-Compilers* will be greyed out.


Click *Install* to continue.

QUARC





A progress bar and embedded video should appear on the installation screen.

QUARC



The QUARC installer automatically configures the Quanser License Manager before completing the installation.




Start MATLAB, and go to *MATLAB Preferences*. Under the *Help* section, in the *Documentation Location* panel, select *Installed locally* (and for some MATLAB versions, click the *Install Documentation* button). Then click *OK*, and restart MATLAB.

STEP 4 Install QUARC on myRIO via NI MAX


Ensure CompactRIO drivers and its accompanied software are installed as outlined in Step 2. Also make sure QUARC is installed as outlined in Step 3 above. Then connect the myRIO to the computer that has QUARC installed via USB and ensure the main power to the myRIO is connected.

A




- Open NI MAX by launching it via the Windows Start Menu.
- Expand the *Remote Systems* tree, and find the *myRIO* device.
- Look at the *System Settings*, and note the *IP Address* of your myRIO. You will need to use it to setup the QUARC test model in Step 5.

B



- Expand the *myRIO* device, and select *Software*.
- Right-click on *Software*, and click on *Add/Remove Software*.


C



Once the *Log In* window pops up, use the login information from your myRIO manual to proceed.


Tip: The default *User name* is *admin*, with no *Password*.

D




When the *LabVIEW Real-Time Software Wizard* dialog pops up, click on *Custom software installation*, and then click *Next*».

E




This warning dialog will appear, click *Yes* to accept the warning.

F




- Select *QUARC for myRIO* from the list.
- From the feature drop-down menu, choose *Install the feature*.
- Click *Next*.

G



A dialog showing the selected items will be displayed.
Click *Next >>* to continue.


H



A progress bar should appear on the installation screen while NI MAX installs QUARC to the myRIO.

Note: The myRIO will reboot a few times during the QUARC installation.

I



When the installation is complete, a message will indicate that the myRIO has been updated successfully. Click *Finish* to exit the dialog window.

My System	
> Data Neighborhood	
> Devices and Interfaces	
> Scales	
> Software	
Remote Systems	
> Ni-myRIO-1900-030546e9	
> Devices and Interfaces	
> Software	
LabVIEW Control Design and Simu	
LabVIEW myRIO Toolkit 16.0.0	
LabVIEW Real-Time 16.0.0	
Network Streams 16.0	
Network Variable Engine 16.0.0	
Ni System Configuration Remote ?	
> NI Web-based Configuration and I	
NI Wireless Certificate Managemei	
> NI-RIO 16.0	
> NI-VISA 16.0.0	
QUARC for myRIO 2.7.2477	
Run-Time Engine for Web Services	

Name	Value
Description	QUARC for myRIO contains the components ...
Version	2.7.2477

- You should see *QUARC for myRIO* in the list of software.
- You can now exit NI MAX.



STEP 5 Software-only Test

The *QUARC Sine and Scope Demo* used in this section is to confirm that QUARC has been installed properly on both the host PC and on the myRIO.

A



Make sure the myRIO is powered (the power LED is lit up).



- Launch MATLAB.
- Open the *QUARC preferences* by typing **quarc_preferences_dialog** in the MATLAB Command Window.
- The *QUARC Preferences* dialog window should appear.
- On the **Model** tab, set the following:
Target type: *linux_rt_armv7*
Default model URI:
tcpip://MYRIO_IP_address:17001

where *MYRIO_IP_address* is the IP address of your myRIO, which can be found via NI MAX (refer to Step 4A).

- On the **Target Type** tab, set the following:
Default target type: *linux_rt_armv7*
Default compiler: *gcc*
- Then press *OK* to close the dialog.

C

QUARC Demo EXAMPLES

QUARC allows you to run Simulink models in real-time. You can create and control the real-time execution entirely through Simulink. QUARC supports multiple targets, including Windows and Linux-based targets.

Quanser products

Basic Features

QUARC Sine and Scope Demo Uses: QUARC, Simulink	Model
QUARC Data Logging Demo Uses: QUARC, Simulink	Model
QUARC Computation Time Demo Uses: QUARC, Simulink	Model
QUARC System Timebase Demo Uses: QUARC, Simulink	Model
QUARC Multi-Rate Demo Uses: QUARC, Simulink	Model

quarc_sine_scope_demo.mtl

Open this model

QUARC Sine and Scope Demo

This is a simple yet important example that can be used to ensure correct operation of QUARC without requiring any special hardware such as a data acquisition card. A sine wave with amplitude 1 and frequency 1 rad/sec is generated and written directly to a scope where it can be traced.


QUARC Sine and Scope Demo

Download here
for instructions

Sine Wave Scope

- Open the QUARC built-in examples by typing **qc_show_demos** in the MATLAB Command Window.
- The *QUARC Demo Examples* Supplemental Software Help window should appear.
- Click on the **QUARC Sine and Scope Demo** under the *Basic Features* category to open the example page.
- On the top-right corner of the *QUARC Sine and Scope Demo* example page, click on **Open this model**.


D



On the taskbar, right click on the QUARC tray icon and select *Target->Remote...*

The *Select remote target* dialog should appear.

E




On the *Select remote target* dialog, make sure to set the Target URI to:

tcpip://MYRIO_IP_address:17000


where *MYRIO_IP_address* is the IP address of your myRIO, which can be found via NI MAX (refer to Step 4A). Then press *OK* to close the dialog. The QUARC Tray icon should now appear green.

F




- Click on the **Monitor & Tune** button on the Hardware tab of the Simulink Toolstrip.
- Double-click on the **Scope** block.
- A sine wave of amplitude 1 should be plotted in real-time.

G



- Double-click on the **Sine Wave** block.
- Change the **Amplitude** to 5, then press the **Apply** button.
- The sine wave's amplitude would be changed to 5 in real-time.
- Change the **Frequency (rad/sec)** to 2, then press the **Apply** button
- The sine wave's frequency would be doubled in real-time

H



Click on the **Stop** button on the Hardware tab of the Simulink Toolstrip to stop the running model.

STEP 6 DAQ Test


The *QUARC Analog Loopback Demo* used in this section is to confirm QUARC has been installed properly. It also tests the myRIO is working properly. Note, it is assumed that the instructions in Step 5 have been followed.

A



Make sure the **myRIO** is powered (the power LED is lit up).

B



Using a jumper wire, connect the **Analog Output Channel #0 (AO0)** to the **Analog Input Channel #0 (AIO)** on the **MXP A connector** on the myRIO.

C

QUARC Demo EXAMPLES

QUARC allows you to run Simulink models in real-time. You can create and control the real-time execution entirely through Simulink. QUARC supports multiple targets, including Windows and Linux-based targets.

Quanser products

Basic Features

QUARC Sine and Scope Demo Uses: QUARC, Simulink	Model
QUARC Data Logging Demo Uses: QUARC, Simulink	Model
QUARC Computation Time Demo Uses: QUARC, Simulink	Model
QUARC System Timebase Demo Uses: QUARC, Simulink	Model
QUARC Multi-Rate Demo Uses: QUARC, Simulink	Model

quarc_analog_loopback_demo.mdl

Open this model

QUARC Analog Loopback Demo


This example is a simple analog loopback test, but it demonstrates a number of important features of QUARC, such as the ability to access and quickly change hardware, multiple targets, online parameter tuning, data streaming, MATLAB logging, using Model Explorer and potentially more.

Download Software
for Instructions




- Open the QUARC built-in examples by typing **qc_show_demos** in the MATLAB Command Window.
- The *QUARC Demo Examples* Supplemental Software Help window should appear.
- Click on the **QUARC Analog Loopback Demo** under the *Using Hardware* category to open the example page.
- On the top-right corner of the *QUARC Analog Loopback Demo* example page, click on **Open this model**.

D




Double-click on the QUARC **HIL Initialize** block.

E




- In the *Board type* option list under the *Main* tab, select **ni_myrio_1900**.
- Click on the **OK** button to close the *HIL Initialize* dialog.

F




- Click on the **Monitor & Tune** button on the Hardware tab of the Simulink Toolstrip.
- Double-click on the **Scope** block.

G




A 1-Hz half-wave rectified sine wave of amplitude 1 should be plotted in real-time. If not, go to the *Troubleshooting* section.

In the model, the Simulink Signal Generator block applies a 1-Hz sine wave signal of a 1-Volt amplitude to the myRIO analog output channel #0. However analog output channel #0 is a 0 to 5V analog output, thus it would not output negative voltage. Due to the jumper wire connection, this signal is then acquired and read through the myRIO analog input channel #0 (note that analog input channel #0 is also a 0 to 5V analog input).



- Double-click on the **Signal Generator** block.
- Change the **Amplitude** to 2, then press the **Apply** button.
- The sine wave's amplitude would be changed to 2 in real-time.
- Change the **Frequency** to 2, then press the **Apply** button
- The sine wave's frequency would be doubled in real-time



Click on the **Stop** button on the Hardware tab of the Simulink Toolstrip to stop the running model.

TROUBLESHOOTING

Review the following recommendations before contacting Quanser's technical support engineers.

The "myRIO Cross-Compilers" is greyed out in features selection during QUARC install.	<ul style="list-style-type: none"> • Ensure your QUARC license file includes the myRIO feature. • Ensure you have installed NI MAX prior to installing QUARC. Note that you can continue installing QUARC without NI MAX, however support for myRIO will not be available. In this case, install NI MAX after installing QUARC, then re-run the QUARC web installer again, and choose Change to allow you to select the myRIO Cross-Compilers feature to be installed.
The myRIO does not appear in the NI Measurement & Automation Explorer (MAX).	<ul style="list-style-type: none"> • Refer to the National Instruments troubleshooting guide for the NI CompactRIO device at www.ni.com/getting-started/set-up-hardware/compactrio/troubleshoot-max.
When running the Analog Loopback Demo, the Scope does not display a sine wave.	<ul style="list-style-type: none"> • Ensure the jumper loopback connection is made on the MXP A connector, as described in Step 6A. • Verify that ni_myrio_1900 was selected in the HIL Initialize dialog, as described in Step 6E.
Getting error: <i>'Error occurred while executing External Mode MEX-file 'quarc_comm': An operating system specific kernel-level... driver for the specified card could not be found. The card or driver may not be installed...'</i> when building a model	<ul style="list-style-type: none"> • Ensure the USB cable is properly connected between the host PC and the myRIO. • Ensure QUARC is installed on the myRIO as outlined in Step 4.

LEARN MORE

To find out about the full range of Quanser control experiments, visit www.quanser.com