



# QBot Platform

## Power User Manual

*v 1.0 – 29th February 2024*

For more information on the solutions Quanser Inc. offers,  
please visit the web site at: <http://www.quanser.com>

Quanser Inc. info@quanser.com  
119 Spy Court Phone : 19059403575  
Markham, Ontario Fax : 19059403576  
L3R 5H6, Canada printed in Markham, Ontario.

This document and the software described in it are provided subject to a license agreement. Neither the software nor this document may be used or copied except as specified under the terms of that license agreement. Quanser Inc. grants the following rights: a) The right to reproduce the work, to incorporate the work into one or more collections, and to reproduce the work as incorporated in the collections, b) to create and reproduce adaptations provided reasonable steps are taken to clearly identify the changes that were made to the original work, c) to distribute and publicly perform the work including as incorporated in collections, and d) to distribute and publicly perform adaptations. The above rights may be exercised in all media and formats whether now known or hereafter devised. These rights are granted subject to and limited by the following restrictions: a) You may not exercise any of the rights granted to You in above in any manner that is primarily intended for or directed toward commercial advantage or private monetary compensation, and b) You must keep intact all copyright notices for the Work and provide the name Quanser Inc. for attribution. These restrictions may not be waived without express prior written permission of Quanser Inc.

#### Waste Electrical and Electronic Equipment (WEEE)



This symbol indicates that waste products must be disposed of separately from municipal household waste, according to Directive 2002/96/EC of the European Parliament and the Council on waste electrical and electronic equipment (WEEE). All products at the end of their life cycle must be sent to a WEEE collection and recycling center. Proper WEEE disposal reduces the environmental impact and the risk to human health due to potentially hazardous substances used in such equipment. Your cooperation in proper WEEE disposal will contribute to the effective usage of natural resources.



**This equipment is designed to be used for educational and research purposes and is not intended for use by the public.** The user is responsible to ensure that the equipment will be used by technically qualified personnel only. While the end-effector board provides connections for external user devices, users are responsible for certifying any modifications or additions they make to the default configuration.

## Table of Contents

A. Charging the battery	3
B. Connecting the battery to the QBot Platform	6
C. Turning ON the QBot Platform	6
D. Turning OFF the QBot Platform	7
E. Low-battery and auto-shutdown	7

## A. Charging the battery



**Caution:** If a battery has been left uncharged for a long time, pay close attention when attempting to charge it again.

1. Power the Optimate battery charger (Figure 1a), that has an SAE charging cable attached.
2. **Skip this step unless your battery becomes disconnected from the battery connectors.** Connect the SAE battery connector cable (highlighted in the red rectangle) to the battery by connecting the red spade connector to the red terminal on the battery and the black spade connector to the black terminal on the battery (Figure 1b).  
**NOTE:** You should only do this step if the connectors become disconnected to the battery. There are other methods of attaching to your battery and charging it that come with this battery charger (Figure 1b). Any of the connector cables can be used, but for day-to-day ease of use we suggest using the SAE cable already connected to the batteries provided.
3. Plug the SAE connector on the cable attached to the battery into the SAE connector on the cable attached to the battery charger.  
**NOTE:** If the battery is connected incorrectly the reverse polarity symbol on the battery charger will light up to let you know (Figure 1c).
4. Charging should start automatically.  
**NOTE:** If the battery is completely drained the internal protection will trigger (Figure 1d). In this case you will have to press the black **Tune** button (Figure 1e) until the red charging button stops flashing and is solid red, this will start the battery charging as the charger will not automatically sense that there is a battery connected.
5. Charging is complete when the LED has moved up to the green check marked section (See Figure 1f). Battery will charge much faster to a 75% charge (< 1 hour). A full charge should take no more than 6 hours.
6. After charging is complete, unplug the SAE connector and re-plug that connector back into the corresponding SAE connector on the QBot Platform. Then disconnect the battery charger from the outlet.

For more information about the battery charger and battery safety please consult the Optimate lithium battery charger manual included with your charger or see the manual [here](#). More information on the Dakota Lithium Iron Phosphate (LFP or LiFePO<sub>4</sub>) batteries can be found [here](#). Some battery parameters have been summarized in Table 1.

#	Item	Value
1	Battery capacity	7 Ah
2	Max continuous discharge	10 A
3	Instantaneous discharge	50 A @ 300ms pulse
4	Connector on battery side	SAE
5	Maximum charge voltage (recommended)	14.4 V
6	Nominal voltage	12.8 – 13.4 V
7	Minimum voltage (cutoff/recommended)	9.0 / 11.0 V
8	Battery weight	907.18 grams
9	Battery dimensions (L x W x H)	15.1 cm x 6.5 cm x 9.5 cm

Table 1. Dakota Lithium 12V 7Ah LiFePO<sub>4</sub> battery characteristics



a. Charger/Balancer



b. Battery & Battery Connector Cables



c. Reverse Polarity Warning



d. Fully Drained Battery or Not Connected



e. Fully Drained Battery/Tune



f. Fully Charged

Figure 1 (continued). Using the batteries and battery charger

## B. Connecting the battery to the QBot Platform

Once the LiFeP04 battery is charged, place it in the QBot Platform by following these steps in Figure 2:

1. Open the wings of the QBot Platform gently tugging on the sides of the wings until the strength of the magnet is released and they pop open.
2. Place the LFP battery in any bay of choice (or both) and connect the SAE connector of the battery cable to the QBot. Guide the wire under the wing rails.
3. Close the wings and you are ready to go.



a. Empty battery bay (right bay #1)



b. battery bay with LFP battery connected to QBot Platform

Figure 2. QBot Platform battery compartment with connectors

## C. Turning ON the QBot Platform

After the batteries are connected, turn on the QBot platform by firmly but gently pressing on the silver power button shown in Figure 3.



Figure 3. Silver power button and emergency motor shutoff

The power LED strips around the QBot Platform should turn red. The LCD will display 'QUANSER INNOVATE EDUCATE' as shown in Figure 4a. After a few seconds, the LCD will

change to display the serial number and say 'qbp-XXXXX starting up...' as shown in Figure 4b. Once the QBot Platform boots up completely, it will switch to a default mode that shows which battery the robot is running on as well as the current battery level, Ethernet/Wi-Fi status as well the IP address that the robot is connected to, an example of which is shown in Figure 4c. It may take several minutes for the robot to connect and display the IP address after the battery level is shown.



a. First bootup message



b. Second bootup message



c. Default display message

Figure 4. LCD startup and default messages

## D. Turning OFF the QBot Platform

To turn OFF the QBot Platform, firmly but gently press on the silver power button once to turn off the robot. **This is the recommended process to initiate shutdown normally.** Please wait for the QBot Platform to fully shut down before disconnecting the battery or tipping the robot on its side (as it will damage the lidar).



Figure 5. Normal shutdown message

If for some reason the QBot Platform will not shut down by clicking the power button, press and hold the power button for approximately 3-4 seconds. The QBot Platform's lights will go off indicating the robot has shut down, however the LCD will hold its last display until you release holding the button. Upon releasing the power button, the power to the QBot Platform computer will be cut off immediately and the LCD will go blank, equivalent to a hard shutdown. This is **not recommended** and should not be required during normal operation, as it may corrupt the SD card on the embedded computer.

## E. Low-battery and auto-shutdown

When the Lithium Iron Phosphate (LFP or LiFePO<sub>4</sub>) batteries get low, safety measures have been built into the QBot Platform to ensure proper shutdown of the robot and maintain the health of the batteries.



### Single battery configuration:

On power up, the QBot Platform will continue using the battery until its voltage goes under 12.0V. At this point, a LOW BATX message will get displayed on the LCD as shown in Figure 6a. Below 11.5V, the QBot Platform will automatically shutdown to prevent battery damage.

### Dual battery configuration:

On power up, the QBot Platform will select the battery with a lower but valid voltage. A valid voltage is any value above 12.0V. When the selected battery voltage goes under 12.0V, the QBot will automatically switch to the other battery so long as it is still valid. This switchover is automatic via an onboard relay and no reboot is required.

If both batteries reach a voltage under 12.0V, the message in Figure 6b will be displayed. At this point, both batteries are used and drained simultaneously until they reach 11.5V to prolong their life. Under 11.5V, the QBot Platform will automatically shut down, displaying the shutdown message. It will first attempt a normal shutdown, but if it is unable to do so, it will disconnect the power after a short period and the LED strips will flash red.



a. Low battery warning message



b. Recharge Message



c. Automatic Shutdown

Figure 6. Low voltage protection and automatic shutdown

© Quanser Inc., All rights reserved.



Solutions for teaching and research. Made in Canada.